

A GUIDE TO THE NEW GENEMARK TISSUE SAMPLING SYSTEM

- HOW TO LOAD THE ALLFLEX TISSUE APPLICATOR
- TAKING A SAMPLE
- RECORDING SAMPLING DATA



Equipment needed

To take a tissue sample using the Allflex Tissue Sampling Unit, the following tools are needed:

1. An Allflex tissue applicator
2. Sample detail card
3. Allflex tissue punches
 - Animals should be restrained during sampling (for example, by packing tightly into a race or using a head bail).
 - Somewhere sheltered to stand for the person writing down animal IDs is also needed.

HOW TO USE THE ALLFLEX TISSUE APPLICATOR

LOADING THE ALLFLEX TISSUE APPLICATOR

1. Remove a punch from the box supplied. Punches should be assembled as one piece; if they have come apart, reassemble.



2. Ensure the tube retainer at the base of the applicator gun is open. Twist the retainer base until the retainer is open.



3. Insert the punch into the tissue applicator as shown below. Twist the tube retainer to lock the punch into the device.



4. Carefully squeeze the applicator handles together, guiding the punch tip into place if necessary. When fully seated the gun bolt will rest flush against the red plastic clip.



5. Release the handles.



6. Remove the red plastic clip by pushing it downward until it unclips. TAKE CARE NOT TO CUT FINGERS ON THE METAL CUTTER – IT IS VERY SHARP.



7. The applicator is now loaded and ready to take a sample. All parts are for one-off use only; do not re-use metal cutters.

TAKING A SAMPLE


1. Check the tissue applicator is loaded with an unused tissue punch (the red plunger is visible in used punches).



Un-used punch

Used punch

2. Slide the gun over the ear, and position the cutter at least 1cm to 2cm from the edge of the ear, taking care to avoid any obvious veins and ridges.



3. Squeeze handles together to take a sample and then release to free the ear. Do this in one swift, fluid, motion.



4. Remove the punch from the device and check that sampling has been successful. If not, discard the punch and re-sample with a new punch.



5. Remove the used cutter from the applicator by pulling the handles apart. This will loosen the cutter. Discard safely.

CAUTION: THE CUTTER IS VERY SHARP, MIND YOUR FINGERS.



RECORDING SAMPLING DATA

1. Be sure to record the following information:

Name and participant (PTPT) code.

2. Sample date. Important because cow IDs may be reused if an animal leaves the herd.

3. Cow numbers/birth IDs. Circle one. Samples cannot be processed unless clearly identified.

4. You will notice each punch has a unique ID and barcode written on the base of the tube. The last five digits are written in larger font on the side of the tube. This is to allow you to match the tissue punch number to the card number. For example, Cow 156 was sampled using punch A 10 02068.

Name	
Participant Code	
Sample Date	
Cow Numbers	Birth IDs (circle)
	A 10 02060
	A 10 02061
	A 10 02062
	A 10 02063
	A 10 02064
	A 10 02065
	A 10 02066
	A 10 02067
156	A 10 02068
	A 10 02069

Full birth IDs should be recorded.

Make sure animals with different years and PTPT codes are clearly noted.

If used, make sure cow numbers are clearly noted.

Make sure sires/males are clearly noted.

As each box of 10 samples is completed and the submission card on the lid is filled out, the lid can be removed by tearing along the perforated hinge line.

Please be sure to send the lid of the box only – the entire box will not fit in the FREEPOST bags provided by GeneMark.

Name <i>A C Farmer</i>	
Participant Code <i>KFWU</i>	
Sample Date <i>6-11-2010</i>	
Cow Numbers	Birth IDs (circle)
	A 10 02060
<i>56</i>	A 10 02061
<i>25</i>	A 10 02062
<i>181</i>	A 10 02063
<i>62</i>	A 10 02064
<i>KFWU-08-21</i>	A 10 02065
<i>PDVD-06-271</i>	A 10 02066
<i>Tag 150</i>	A 10 02067
<i>KFWU-06-22 (Bull)</i>	A 10 02068
	A 10 02069

- Record sample information on the lid of the sample boxes (as shown to the right of this page).
- Be sure to assign the correct animals to the correct sample number.
- After sampling is complete, samples should be placed into a ziplock bag and into the freezer until you are ready to courier them to GeneMark (140 Riverlea Rd, Hamilton).
- Note: The samples are biological biopsies. They will survive several days in the post but should be treated with care prior to postage.
- DURING SAMPLING KEEP THE SAMPLES OUT OF THE SUN.
- AFTER SAMPLING FREEZE THEM UNTIL READY TO COURIER TO GENEMARK. FAILURE TO DO THIS MAY CAUSE ROTTING AND TEST FAILURES.
- DO NOT STORE SAMPLES AT ROOM TEMPERATURE.

The instructions on these pages can be viewed on video. Contact your Customer Relationship Manager or watch it on the LIC website: <http://www.lic.co.nz/genemark>.